# Bingzun Infrared Thermometer

**Brand Name:** Bing Zun

**Model:** BZ-R6

**Features:**

- FDA Certificate
- Class I Medical Device(lowest risk)
- General Hospital Use
- 510K Exempt
- CE Certificate
- ISO Certificate
- FCC Certificate

# Product Pictures

Henan Bingzun Industrial Co., Ltd  
**Bing Zun**

**1** second  
Temperature measurement

**Infrared Thermometer R6**

Non-contact Infrared Thermometer  
TEMPERATURE | INSTANT | ACCURATE

**Bing Zun Infrared  
Thermometer  
User Manual**

The manufacturer has reservation right of  
product technical index without prior notice.  
The manufacturer has the right to modify the  
technical information of the product without  
prior notice.

**Product certificate**

Product name: Infrared thermometer

Model number: R6

Product grade: Qualified product

Executive standard: GB/T 2015-2458

Production date:

Henan Bingzun Industrial Co., Ltd

Address: 10 meters west to the south of  
Gufu Road, Xin'an Industry Cluster District,  
Luoyang City, Henan Province, China

# Product Packaging

BOX\*1  
THERMOMETER\*1  
USER MANUAL\*1  
PRODUCT CERTIFICATE\*1

# Master Case Packaging

40PC/CTN

CARTON SIZE:

- 49.5 \* 40.6 \* 29 CM
- 19.5 \* 16.0 \* 11.5 INCH

NET WEIGHT: 6.12KG/13.5LB/CTN

GROSS WEIGHT: 7.02KG/15.5LB/CTN

# Product Spec Sheet

BINGZUN BZR6


Delicate, efficient and convenient

Product name	Medical infrared thermometer	
Product model	E1501	
Screen specification	4 * 14 (TN full transparent positive display, 28.5mm * 34.3mm * 2.8mm)	
Probe type	Tts-m5	
product type	Human body temperature measuring device	
Product weight	100g (without battery)	
Measuring distance	3CM-5CM	
Measurement accuracy	±0.2° C(0.5° F)	
Automatic shut off	About 20s	
Measuring time	0.5s	
Supply voltage	DC 3V (2 No.7 alkaline batteries / in series)	
Frequency of use	Press more than 10000 times	
Show exact digits	0.1° C (0.1° F)	
Storage temperature	-20-55 degrees C	
Operating temperature	5 ° C ~ 40 ° C, optimal temperature 25 °	
Relative humidity	≤85%	
Temperature range	32.0-42.99C	
High measurement accuracy (using plug-in only ADC algorithm processing)	32~35.9° C (93.2~96.6° F)	±0.3° C(0.5° F)
	36~39° C (96.8~102.2° F)	±0.2° C(0.4° F)
	39~43° C (102.2~109.4° F)	±0.3° C(0.5° F)

**Product structure:**

**Delicate, efficient and convenient**

1.	Infrared sensor
2.	Liquid crystal (LCD) display
3.	Backlight on / off button
4.	Up button (view historical data)
5.	Down button (view historical data)
6.	Sound switch button
7.	Measuring switch
8.	Handle
9.	Battery cover

**Reference temperature:**

Measuring location	Normal temperature (° C)	Normal Fahrenheit (° f)
Anus	36.6~38	97.8~100.4
oral cavity	35.5~37.5	95.9~99.5
Armpit	34.7~37.3	94.4~99.1
Ear	35.8~38	96.4~100.4
forehead	36~37.2	97.4~98.4

# Product Certificates

BINGZUN BZR6

APPROVED

- Home
- Food
- Drugs
- Medical Devices
- Radiation-Emitting Products
- Vaccines, Blood & Biologics
- Animal & Veterinary
- Cosmetics
- Toba

# Establishment Registration & Device Listing

FDA Home Medical Devices Databases

<a href="#">New Search</a>	<a href="#">Back To Search Results</a>
<b>Proprietary Name:</b>	Non-contact Infrared Thermometer; BZ-R3, R5, R6, R7, R8, R9, R11
<b>Classification Name:</b>	THERMOMETER, CLINICAL COLOR CHANGE
<b>Product Code:</b>	<u>FQZ</u>
<b>Device Class:</b>	<b>1</b>
<b>Regulation Number:</b>	<u>880.2900</u>
<b>Medical Specialty:</b>	<b>General Hospital</b>
<b>Registered Establishment Name:</b>	<u>HENAN BINGZUN INDUSTRIAL CO., LTD</u>
<b>Owner/Operator:</b>	<u>Henan Bingzun Industrial Co., Ltd</u>
<b>Owner/Operator Number:</b>	10073151
<b>Establishment Operations:</b>	Manufacturer

FDA Class I medical device:  
Device with lowest risk  
to the patient and/or user  
Pregnancy and Baby Safe

# Product Classification

[FDA Home](#) [Medical Devices](#) [Databases](#)

[New Search](#)

[Back to Search Results](#)

<b>Device</b>	Thermometer, Clinical Color Change
<b>Regulation Description</b>	Clinical color change thermometer.
<b>Regulation Medical Specialty</b>	<b>General Hospital</b>
<b>Review Panel</b>	General Hospital
<b>Product Code</b>	<b>FQZ</b>
<b>Premarket Review</b>	<a href="#">Gastrorenal, ObGyn, General Hospital, and Urology Devices (OHT3)</a> <a href="#">Drug Delivery and General Hospital Devices, and Human Factors (DHT3C)</a>
<b>Submission Type</b>	<b>510(K) Exempt</b> 
<b>Regulation Number</b>	<a href="#">880.2900</a>
<b>Device Class</b>	<b>1</b>
<b>Total Product Life Cycle (TPLC)</b>	<a href="#">TPLC Product Code Report</a>
<b>GMP Exempt?</b>	No
<b>Summary Malfunction Reporting</b>	Eligible

**510(K) Exempt:  
No 510K Needed  
No Pre-market Review Need**

**Note:** FDA has exempted almost all class I devices (with the exception of reserved devices) from the premarket notification requirement, including those devices that were exempted by final regulation published in the *Federal Registers* of December 7, 1994, and January 16, 1996. It is important to confirm the exempt status and any limitations that apply with [21 CFR Parts 862-892](#). Limitations of device exemptions are covered under 21 CFR XXX.9, where XXX refers to Parts 862-892.

If a manufacturer's device falls into a generic category of exempted class I devices as defined in [21 CFR Parts 862-892](#), a premarket notification application and fda clearance is not required before marketing the device in the U.S. however, these manufacturers are required to register their establishment. Please see the [Device Registration and Listing website](#) for additional information.

### Recognized Consensus Standard

- 6-70 ASTM E825-98 (Reapproved 2016)  
[Standard Specification for Phase Change-Type Disposable Fever Thermometer for Intermittent Determination of Human Temperature](#)

<b>Implanted Device?</b>	No
<b>Life-Sustain/Support Device?</b>	No
<b>Third Party Review</b>	Not Third Party Eligible

# CE CERTIFICATE

## Certificate Of Compliance

Certificate Number: UNIA20040907EC-01

**Applicant** : Henan BingZun Industrial Co., Ltd  
50 meters next to the south of Guihua Villa, Xin'an Industry Cluster District, Luoyang City, Henan Province, China

**Manufacturer** : Henan BingZun Industrial Co., Ltd  
50 meters next to the south of Guihua Villa, Xin'an Industry Cluster District, Luoyang City, Henan Province, China

**Product** : Non-contact Electronic Thermometer

**Trade Name** : Bingzun, Cry, Yostand, Rainbow+

**M/N** : BZ-R6, R3, R5, R6, R7, R8, R9, R11, R12, R13, ET01, ET03, ET05, TWQ-01, Z3, Z5, Z6

**Test Standard:** EN 61326-1:2013

The EUT described above has been tested by us with the listed standards and found in compliance with the council EMC Directive 2014/30/EU. It is possible to use CE marking to demonstrate the compliance with EMC Directive.

The certificate applies to the tested sample above mentioned only and shall not imply an assessment of the whole production. It is only valid in connection with the test report number: UNIA20040907ER-01.

2020-04-22

Date of issue

Hoffer Lau

# CE

**Shenzhen United Testing Technology Co., Ltd.**

2F, Annex Bldg, Jiahuangyuan Tech Park, #365 Baotian 1 Rd, Tiegang Community, Xixiang Str, Bao'an District, Shenzhen, China

Tel: +86-755-86180996 Fax: +86-755-86180156

<http://www.uni-lab.hk> [E-mail: hofferlau@uni-lab.hk](mailto:hofferlau@uni-lab.hk)

# FCC CERTIFICATE

## SUPPLIER'S DECLARATION OF CONFORMITY

Certificate Number: UNIA20040905EC-01

We Hereby, Declare that the essential requirements set out in the 47 CFR Sections 2.906 of FCC, have been fully fulfilled on our product with indication below:

### Applicant's info

Company : Henan BingZun Industrial Co., Ltd  
Address : 50 meters next to the south of Guihua Villa, Xin'an Industry Cluster District, Luoyang City, Henan Province, China

### Product Info

Product Name : Non-contact Electronic Thermometer  
Model No. : BZ-R6, R3, R5, R6, R7, R8, R9, R11, R12, R13, ET01, ET03, ET05, TWQ-01, Z3, Z5, Z6  
Trade Mark : Bingzun, Cry, Yostand, Rainbow+

### Applied Standard

Standards  
E-CFR Title 47 Part 15 Subpart B

### Responsibilities and obligations

Strictly follow the FCC rules of 2.906, 2.908, 2.909.

### Responsible party in US.

Company: \_\_\_\_\_ Email: \_\_\_\_\_ Signature: \_\_\_\_\_  
Address: \_\_\_\_\_ Tel: \_\_\_\_\_ Date: \_\_\_\_\_  
Contact person (print name): \_\_\_\_\_ Fax: \_\_\_\_\_

### Test Firm Used to Take Measurements

Name: Shenzhen United Testing Technology Co., Ltd.  
Address: #365 Baotian 1 Rd, Xixiang Str, Bao'an District, Shenzhen, Guangdong, China  
Contact Name: Mr Liu ze  
Tel: +86-755-86180996  
Fax: +86-755-86180156  
Email: hofferlau@uni-lab.hk  
Designation No.: CN1227  
FRN: 0027159896

# ISO13485:2016

Medical Device Management

CERTIFICATION  
EUROPE HONG KONG

## Certificate of Registration

Certificate NO: 10136589M

This is to Certify that the Medical Devices Industry Management System of

**Henan Bingzun Industrial Co., Ltd.**

50 Meters next to the south of Guihua Villa, Xin'an Industry  
Cluster District, Luoyang City, Henan Province, China

Has been audited to the following Medical Devices Industry Management System standard:

**ISO13485:2016**

This system is valid for the

**Production and sales of non-contact electronic thermometer, digital blood pressure monitor, ear scope, electronic thermometer, atomizer, medical flashlight, ear scope, pulse oximeter, sputum aspirator, blood pressure meter, tourniquet, diagnostic hammer, stethoscope, skin nerve test needle, mercury-free thermometer**

Date of issue: Mar. 17, 2020

Date of expiry: Mar. 16, 2023

This certificate will not remain valid only if the certified organization accepts at least one surveillance audit annually within the validity period of the certificate in which the surveillance audit conforming mark is in the designated position on the certificate.

Issued by:

*Franco*

12months

24months

CERTIFICATION EUROPE(HONGKONG)LIMITED WWW.CCE-HK.COM

# Test Report

- **EMC TEST REPORT**
- **CHINA MEDICAL DEVICE TEST REPORT**

# EMC TEST REPORT

**Product:** Non-contact Electronic Thermometer  
**Trade Name:** Bingzun, Cry, Yostand, Rainbow+  
**Model Name:** BZ-R6  
**Serial Model:** R3, R5, R6, R7, R8, R9, R11, R12, R13,  
ET01, ET03, ET05, TWQ-01, Z3, Z5, Z6  
**Report No.:** UNIA20040907ER-01

## Prepared for

Henan Bingzun, Cry, Yostand, Rainbow+ Industrial Co., Ltd  
50 meters next to the south of Guihua Villa, Xin'an Industry Cluster District,  
Luoyang City, Henan Province, China

## Prepared by

Shenzhen United Testing Technology Co., Ltd.  
2F, Annex Bldg, Jiahuangyuan Tech Park, #365 Baotian 1 Rd, Tiegang  
Community, Xixiang Str, Bao'an District, Shenzhen, China

## TEST RESULT CERTIFICATION

**Applicant's name**.....: Henan Bingzun, Cry, Yostand, Rainbow+ Industrial Co., Ltd  
**Address**.....: 50 meters next to the south of Guihua Villa, Xin'an Industry Cluster District, Luoyang City, Henan Province, China

**Manufacture's Name**.....: Henan Bingzun, Cry, Yostand, Rainbow+ Industrial Co., Ltd  
**Address**.....: 50 meters next to the south of Guihua Villa, Xin'an Industry Cluster District, Luoyang City, Henan Province, China

**Product description**

**Product name**.....: Non-contact Electronic Thermometer  
**Trade Mark**.....: Bingzun, Cry, Yostand, Rainbow+  
**Model and/or type reference** ..: BZ-R6, R3, R5, R6, R7, R8, R9, R11, R12, R13, ET01, ET03, ET05, TWQ-01, Z3, Z5, Z6

**Standards**.....: EN 61326-1:2013

This device described above has been tested by Shenzhen United Testing Technology Co., Ltd., and the test results show that the equipment under test (EUT) is in compliance with the EMC Directive 2014/30/EU requirements. And it is applicable only to the tested sample identified in the report.

This report shall not be reproduced except in full, without the written approval of UNI, this document may be altered or revised by Shenzhen United Testing Technology Co., Ltd., personnel only, and shall be noted in the revision of the document.

**Date of Test**.....:   
**Date (s) of performance of tests**.....: Mar. 13, 2020 ~ Mar. 22, 2020  
**Date of Issue**.....: Mar. 22, 2020  
**Test Result**.....: Pass

Prepared by:

*Bob Liao*

Bob Liao/Editor

Reviewer:

*Kahn Yang*  
 Kahn Yang/Supervisor

Approved & Authorized Signer:

*Liuze*  
 Liuze/Manager

Table of Contents

Page

1 TEST SUMMARY.....	4
2 GENERAL INFORMATION.....	7
2.1 GENERAL DESCRIPTION OF EUT.....	7
2.2 DESCRIPTION OF THE TEST MODES.....	8
2.3 DESCRIPTION OF TEST SETUP.....	8
2.4 DESCRIPTION TEST PERIPHERAL AND EUT PERIPHERAL.....	9
3 RADIATED EMISSIONS MEASUREMENT.....	12
3.1 RADIATED EMISSION LIMIT.....	12
3.2 TEST SETUP.....	12
3.3 TEST PROCEDURE.....	13
3.4 TEST RESULT.....	13
4 EMC IMMUNITY TEST.....	16
4.1 STANDARD COMPLIANCE/SERVIRITY LEVEL/CRITERIA.....	16
4.2 GENERAL PERFORMANCE CRITERIA.....	16
5 ELECTROSTATIC DISCHARGE IMMUNITY TEST (ESD).....	17
5.1 TEST SPECIFICATION.....	17
5.2 TEST SETUP.....	17
5.3 TEST PROCEDURE.....	18
5.4 TEST RESULT.....	19
6 RADIATED, RADIO-FREQUENCY, ELECTROMAGNETIC FIELD IMMUNITY TEST (RS).....	20
6.1 TEST SPECIFICATION.....	20
6.2 TEST SETUP.....	20
6.3 TEST PROCEDURE.....	21
6.4 TEST RESULT.....	22
7 PHOTO OF EUT.....	23
8 PHOTO OF TEST.....	29

# 1 TEST SUMMARY

## TEST RESULTS

Test procedures according to the technical standards:

EMC Emission				
Standard	Test Item	Limit	Judgment	Remark
EN 61326-1:2013	Conducted Emissions From The AC Mains Power Ports	Class B	N/A	
	Radiated Emissions	Class B	PASS	NOTE (1)
EN 61000-3-2:2014	Harmonic Current Emission	---	N/A	NOTE (2)
EN 61000-3-3:2013	Voltage Fluctuations & Flicker	---	N/A	
EMC Immunity				
Section EN 61326-1:2013	Test Item	Performance Criteria	Judgment	Remark
EN 61000-4-2:2009	Electrostatic Discharge	B	PASS	
EN 61000-4-3:2006 +A1:2008+A2:2010	RF Electromagnetic Field	A	PASS	
EN 61000-4-4:2012	Fast Transients	B	N/A	
EN 61000-4-5:2014 /A1:2017	Surges	B	N/A	
EN 61000-4-6:2014 /AC:2015	Injected Current	A	N/A	
EN 61000-4-8:2010	Power Frequency Magnetic Field	A	N/A	
EN 61000-4-11:2004 /A1:2017	Volt. Interruptions Volt. Dips	B / B / C / C	N/A	NOTE (3)

Note:

- (1) If the highest frequency of the internal sources of the EUT is less than 108 MHz, the measurement shall only be made up to 1 GHz.  
 If the highest frequency of the internal sources of the EUT is between 108 MHz and 500 MHz, the measurement shall only be made up to 2 GHz.  
 If the highest frequency of the internal sources of the EUT is between 500 MHz and 1GHz, the measurement shall only be made up to 5 GHz.  
 If the highest frequency of the internal sources of the EUT is above 1 GHz, the measurement shall be made up to 5 times of the highest frequency or 6 GHz, whichever is less.
- (2) The power consumption of EUT is less than 75W and no Limits apply.
- (3) Voltage Dip: 100% reduction – Performance Criteria B  
 Voltage Dip: 100% reduction – Performance Criteria B  
 Voltage Dip: 30% reduction – Performance Criteria C  
 Voltage Interruption: 100% Interruption – Performance Criteria C
- (4) For client's request and manual description, the test will not be executed.
- (5) "N/A" denotes test is not applicable in this Test Report.

## TEST FACTORY

Test Firm : Shenzhen United Testing Technology Co., Ltd.  
Address : 2F, Annex Bldg, Jiahuangyuan Tech Park, #365 Baotian 1 Rd, Tiegang Community, Xixiang Str, Bao'an District, Shenzhen, China

The testing quality ability of our laboratory meet with "Quality Law of People's Republic of China" Clause 19. The testing quality system of our laboratory meets with ISO/IEC-17025 requirements. This approval result is accepted by MRA of APLAC.

Our test facility is recognized, certified, or accredited by the following organizations:

A2LA Certificate Number: 4747.01

The EMC Laboratory has been accredited by A2LA, and in compliance with ISO/IEC 17025:2017 General Requirements for testing Laboratories.

FCC Registration Number: 674885

The EMC Laboratory has been registered and fully described in a report filed with the (FCC) Federal Communications commission.

IC Registration Number: 21947

The EMC Laboratory has been registered and fully described in a report filed with the (IC) Industry Canada.

**MEASUREMENT UNCERTAINTY**

The reported uncertainty of measurement  $y \pm U$ , where expanded uncertainty  $U$  is based on a standard uncertainty multiplied by a coverage factor of  $k=2$ , providing a level of confidence of approximately 95 %.

**A. Conducted Measurement:**

Test Site	Method	Measurement Frequency Range	U, (dB)	NOTE
UNI01	ANSI	9KHz ~ 150KHz	3.18	
		150 KHz ~ 30MHz	2.70	

**B. Radiated Measurement:**

Test Site	Method	Measurement Frequency Range	U, (dB)	NOTE
UNI02	ANSI	9KHz ~ 30MHz	2.50	
		30MHz ~ 200MHz	3.43	
		200MHz ~ 1000MHz	3.57	
		1GHz ~ 6 GHz	4.13	

## 2 GENERAL INFORMATION

### 2.1 GENERAL DESCRIPTION OF EUT

Equipment:	Non-contact Electronic Thermometer	
Trade Mark:	Bingzun, Cry, Yostand, Rainbow+	
Model Name:	BZ-R6	
Serial No.:	R3, R5, R6, R7, R8, R9, R11, R12, R13, ET01, ET03, ET05, TWQ-01, Z3, Z5, Z6	
Model Difference:	All model's the function, software and electric circuit are the same, only with a product color and model named different. Test sample model: BZ-R6.	
Product Description:	The EUT is a Non-contact Electronic Thermometer.	
	Operating frequency:	N/A
	Connecting I/O port:	N/A
	Based on the application, features, or specification exhibited in User's Manual, the EUT is considered as an ITE/Computing Device. More details of EUT technical specification, please refer to the User's Manual.	
Power Source:	DC 3.0V	

## 2.2 DESCRIPTION OF THE TEST MODES

To investigate the maximum EMI emission characteristics generates from EUT, the test system was pre-scanning tested base on the consideration of following EUT operation mode or test configuration mode which possible have effect on EMI emission level. Each of these EUT operation mode(s) or test configuration mode(s) mentioned above was evaluated respectively.

Pretest Mode	Description
Mode 1	Running

For Conducted Test	
Pretest Mode	Description
Mode 1	Running


For Radiated Test	
Pretest Mode	Description
Mode 1	Running

For EMS Test	
Pretest Mode	Description
Mode 1	Running

Note: The test modes were carried out for all operation modes (include link and idle).

## 2.3 DESCRIPTION OF TEST SETUP

Note: The EUT tested system was configured as upper figure, unless otherwise a special operating condition is specified in the following during the testing.

2.4 DESCRIPTION TEST PERIPHERAL AND EUT PERIPHERAL

The EUT has been tested as an independent unit together with other necessary accessories or support units. The following support units or accessories were used to form a representative test configuration during the tests.

Item	Equipment	Mfr/Brand	Model/Type No.	Note
E-1	Non-contact Electronic Thermometer	Bingzun, Cry, Yostand, Rainbow+	BZ-R6	EUT

Item	Shielded Type	Ferrite Core	Length	Note

Note:

1. The support equipment was authorized by Declaration of Confirmation.
2. For detachable type I/O cable should be specified the length in cm in 『Length』 column.
3. “YES” is means “shielded” “with core”; “NO” is means “unshielded” “without core”.

2.5 MEASUREMENT INSTRUMENTS LIST

Item	Equipment	Manufacturer	Model No.	Serial No.	Calibrated until
Conduction Emissions Measurement					
1	Conducted Emission Test Software	EZ-EMC	Ver.CCS-3A1-CE	N/A	N/A
2	AMN	Schwarzbeck	NNLK8121	8121370	2020.10.15
3	AMN	ETS	3810/2	00020199	2020.10.15
4	AAN	TESEQ	T8-Cat6	38888	2020.10.15
5	Pulse Limiter	CYBRTEK	EM5010	E115010056	2020.05.26
6	EMI Test Receiver	Rohde&Schwarz	ESCI	101210	2020.10.15
Radiated Emissions Measurement					
1	Radiated Emission Test Software	EZ-EMC	Ver.CCS-03A1	N/A	N/A
2	Horn Antenna	Sunol	DRH-118	A101415	2020.10.18
3	Broadband Hybrid Antenna	Sunol	JB1	A090215	2020.11.15
4	PREAMP	HP	8449B	3008A00160	2020.10.21
5	PREAMP	HP	8447D	2944A07999	2020.05.26
6	EMI Test Receiver	Rohde&Schwarz	ESBZ-R6	101891	2020.10.15
7	MXA Signal Analyzer	Keysight	N9020A	MY51110104	2020.10.15
8	Active Loop Antenna	Com-Power	AL-310R	10160009	2020.05.28
9	Horn Antenna	Schwarzbeck	BBHA9120D	9120D-1680	2020.05.28
10	Horn Antenna	A-INFOMW	LB-180400-KF	J211060660	2020.10.23
11	Loop Antenna	Beijing daze Technology	ZN30401	13015	2020.10.15
12	EM Clamp	Schwarzbeck	MDS21	03350	2020.10.20
Harmonic / Flicker Measurement					
1	Power Analyzer	California Instrumnets	PACS-1	X71719	2020.10.15
2	AC Power Source	California Instrumnets	5001ix	HK53570	2020.10.15
Electrostatic Discharge Test					
1	ESD Generator	EVERFINE	EMS61000-2A	P185811CA837112 1	2020.10.17
RS Test					
1	Power Meter	Agilent	E4419B	QB4331226	2020.10.10
2	Power Sensor	Agilent	8481A	MY41092622	2020.10.10
3	Power Sensor	Agilent	8481A	US37296783	2020.10.10
4	Signal Generator	Agilent	N5182A	MY46240556	2020.10.10
5	Power Amplifier	MICOTOP	MPA-80-1000-250	1711489	2020.10.10
6	Power Amplifier	MICOTOP	MPA-1000-3000-7 5	1711488	2020.10.10
7	Power Amplifier	MICOTOP	MPA-3000-6000-5 0	MPA1706275	2020.10.10
8	Bilog Antenna	TESEQ	CBL6111D	34678	2020.10.10
9	Horn Antenna	Schwarzbeck	BBHA9120D	9120D-1680	2020.05.28

Item	Equipment	Manufacturer	Model No.	Serial No.	Calibrated until
<b>Electrical Fast Transient / Burst Immunity Test</b>					
1	EMS Test Control System	Shanghai Lioncel	SCU-614AS	SCU614S0160601	N/A
2	EFT/B Generator	Shanghai Lioncel	EFT-404S	EFT404S0160601	2020.10.15
<b>Surge Test</b>					
1	EMS Test Control System	Shanghai Lioncel	SCU-614AS	SCU614S0160601	N/A
2	Surge Generator	Shanghai Lioncel	LSG-506S	LSG506S0160601	2020.10.15
3	CDN	Shanghai Lioncel	CDN-532S	CDN532S0160601	2020.10.15
<b>CS Test</b>					
1	CS	SCHLODER	CDG-6000-25	126A1280/2014	2020.10.10
2	CDN	SCHLODER	CDN-M2+3	A2210275/2014	2020.10.10
3	EM Clamp	SCHLODER	EMCL-20	132A1283	2020.10.10
4	Attenuator	Nemtest	ATT-6DB-100	A100W224	2020.10.10
5	Audio Analyzer	R&S	UPL	100419	2020.10.10
6	Universal Radio Communication Tester	R&S	CMW500	117239	2020.10.10
7	Universal Radio Communication Tester	R&S	CMU200	111764	2020.10.10
8	Audio Analyzer	R&S	UPL	100689	2020.10.10
9	Audio Breakthrough Shielding Box	SKET	SB_ABT/C35	N/A	2020.10.10
10	Ear Simulator	SKET	AE_ABT/C35	N/A	2020.10.10
11	Mouth Simulator	SKET	AM_ABT/C35	N/A	2020.10.10
12	1KHz Standard Source	SKET	MSC_ABT/C35	N/A	2020.10.10
<b>Power-frequency magnetic fields Test</b>					
1	Magnetic Field Test System	Shanghai Lioncel	PMF801C-T	PMF801C-T016070 1	2020.05.26
<b>Voltage dips and interruptions Test</b>					
1	Voltage SAG Simulator	Shanghai Lioncel	VDS-1101	VDS11010160601	2020.10.15
2	Adjustable Power Supply	Shanghai Lioncel	RGL-210	RGL2100151001	N/A

### 3 RADIATED EMISSIONS MEASUREMENT

#### 3.1 RADIATED EMISSION LIMIT

Below 1000MHz:

Frequency (MHz)	Class A		Class B	
	10m	3m	10m	3m
	dBuV/m		dBuV/m	
30~230	40	50	30	40
230~1000	47	57	37	47

Above 1000MHz:


Frequency (MHz)	Class A		Class B	
	PK	AV	PK	AV
	dBuV/m		dBuV/m	
1000~3000	76	56	70	50
3000~6000	80	60	74	54

Note:


1. The tighter limit applies at the band edges.
2. Emission level (dBuV/m)=20log Emission level (uV/m).

#### 3.2 TEST SETUP

##### 1. Radiated Emission Test-Up Frequency Below 1000MHz

## 2. Radiated Emission Test-Up Frequency Above 1000MHz


### 3.3 TEST PROCEDURE

1. The EUT was placed on the top of a rotating table 0.8 meters above the ground at a 3 meter semi-anechoic chamber room. The table was rotated 360 degrees to determine the position of the highest radiation.
2. The EUT was set 3 meters away from the interference-receiving antenna, which was mounted on the top of a variable-height antenna tower.
3. The height of antenna is varied from 1 meter to 4 meters above the ground to determine the maximum value of the field strength. Both horizontal and vertical polarizations of the antenna are set to make the measurement.
4. For each suspected emission, the EUT was arranged to its worst case and then the antenna was tuned to heights from 1 meter to 4 meters and the rotatable table was turned from 0 degrees to 360 degrees to find the maximum reading.
5. The test-receiver system was set to quasi-peak detect function and specified bandwidth with maximum hold mode when the test frequency is below 1GHz.
6. For the actual test configuration, please refer to the related Item EUT Test Photos.

### 3.4 TEST RESULT

PASS


Temperature:	24°C	Relative Humidity:	48%
Test Voltage:	DC 3.0V	Pressure:	1010hPa
Test Mode:	Running	Polarization:	Horizontal

No.	Frequency (MHz)	Reading (dBuV)	Correction factor(dB/m)	Result (dBuV/m)	Limit (dBuV/m)	Margin (dB)	Degree (deg.)	Height (cm)	Remark
1*	30.2111	30.35	-7.31	23.04	40.00	-16.96			peak
2	54.4516	34.75	-21.12	13.63	40.00	-26.37			peak
3	124.5690	26.83	-16.52	10.31	40.00	-29.69			peak
4	452.7197	27.64	-12.26	15.38	47.00	-31.62			peak
5	607.7867	26.98	-9.43	17.55	47.00	-29.45			peak
6	955.4381	34.67	-6.10	28.57	47.00	-18.43			peak

Remark: Absolute Level = Reading Level + Factor, Margin = Absolute Level – Limit  
 Factor = Ant. Factor + Cable Loss – Pre-amplifier

Temperature:	24°C	Relative Humidity:	48%
Test Voltage:	DC 3.0V	Pressure:	1010hPa
Test Mode:	Running	Polarization:	Vertical

No.	Frequency (MHz)	Reading (dBuV)	Correction factor(dB/m)	Result (dBuV/m)	Limit (dBuV/m)	Margin (dB)	Degree (deg.)	Height (cm)	Remark
1*	30.1054	27.11	-7.23	19.88	40.00	-20.12			peak
2	33.4449	28.57	-9.85	18.72	40.00	-21.28			peak
3	54.8348	33.83	-21.16	12.67	40.00	-27.33			peak
4	360.4477	34.97	-14.05	20.92	47.00	-26.08			peak
5	618.5369	27.25	-9.56	17.69	47.00	-29.31			peak
6	955.4381	29.54	-6.10	23.44	47.00	-23.56			peak

Remark: Absolute Level = Reading Level + Factor, Margin = Absolute Level – Limit  
 Factor = Ant. Factor + Cable Loss – Pre-amplifier

## 4 EMC IMMUNITY TEST

### 4.1 STANDARD COMPLIANCE/SERVIRITY LEVEL/CRITERIA

Tests Standard No.	TEST SPECIFICATION	Test Mode Test Ports	Perform Criteria
ESD IEC/EN 61000-4-2	8kV air discharge 4kV contact discharge	Direct Mode	B
	4kV HCP discharge 4kV VCP discharge	Indirect Mode	B
RS IEC/EN 61000-4-3	80 MHz to 1000 MHz, 1400 MHz to 2000 MHz, 2000 MHz to 2700 MHz, 1000Hz, 80%, AM modulated	Enclosure	A

### 4.2 GENERAL PERFORMANCE CRITERIA

According to EN 61326-1 standard, the general performance criteria as following:

Criterion A	The equipment shall continue to operate as intended without operator intervention. No degradation of performance, loss of function or change of operating state is allowed below a performance level specified by the manufacturer when the equipment is used as intended. The performance level may be replaced by a permissible loss of performance. If the minimum performance level or the permissible performance loss is not specified by the manufacturer, then either of these may be derived from the product description and documentation, and by what the user may reasonably expect from the equipment if used as intended.
Criterion B	During the application of the disturbance, degradation of performance is allowed. However, no unintended change of actual operating state or stored data is allowed to persist after the test. After the test, the equipment shall continue to operate as intended without operator intervention; no degradation of performance or loss of function is allowed, below a performance level specified by the manufacturer, when the equipment is used as intended. The performance level may be replaced by a permissible loss of performance. If the minimum performance level (or the permissible performance loss), or recovery time, is not specified by the manufacturer, then either of these may be derived from the product description and documentation, and by what the user may reasonably expect from the equipment if used as intended.
Criterion C	Loss of function is allowed, provided the function is self-recoverable, or can be restored by the operation of the controls by the user in accordance with the manufacturer's instructions. A reboot or re-start operation is allowed. Information stored in non-volatile memory, or protected by a battery backup, shall not be lost.


## 5 ELECTROSTATIC DISCHARGE IMMUNITY TEST (ESD)

### 5.1 TEST SPECIFICATION

Basic Standard:	IEC/EN 61000-4-2
Discharge Impedance:	330 ohm / 150 pF
Required Performance:	B
Discharge Voltage:	Air Discharge: 2kV/4kV/8kV (Direct) Contact Discharge: 2kV/4kV (Direct/Indirect)
Polarity:	Positive & Negative
Number of Discharge:	Air Discharge: min. 20 times at each test point Contact Discharge: min. 200 times in total
Discharge Mode:	Single Discharge
Discharge Period:	1 second minimum

### 5.2 TEST SETUP

Note:

### TABLE-TOP EQUIPMENT

The configuration consisted of a wooden table 0.8 meters high standing on the Ground Reference Plane. The GRP consisted of a sheet of aluminum at least 0.25mm thick. A Horizontal Coupling Plane (1.6m x 0.8m) was placed on the table and attached to the GRP by means of a cable with 940k total impedance. The equipment under test, was installed in a representative system as described in section 7 of IEC /EN 61000-4-2, and its cables were placed on the HCP and isolated by an insulating support of 0.5mm thickness. A distance of 0.8-meter minimum was provided between the EUT and the walls of the laboratory and any other metallic structure.

### FLOOR-STANDING EQUIPMENT

The equipment under test was installed in a representative system as described in section 7 of IEC/EN 61000-4-2, and its cables were isolated from the Ground Reference Plane by an insulating support of 0.1 meter thickness. The GRP was consisted of a sheet of aluminum that is at least 0.25mm thick, and extended at least 0.5 meters from the EUT on all sides.

## 5.3 TEST PROCEDURE

The test generator necessary to perform direct and indirect application of discharges to the EUT in the following manners:

1. Electrostatic discharges were applied only to those points and surfaces of the EUT that are accessible to users during normal operation. The test was performed with at least ten single discharges on the pre-selected points in the most sensitive polarity.

The time interval between two successive single discharges was at least 1 second.

The ESD generator was held perpendicularly to the surface to which the discharge was applied and the return cable was at least 0.2 meters from the EUT.

Contact discharges were applied to the non-insulating coating, with the pointed tip of the generator penetrating the coating and contacting the conducting substrate.

Air discharges were applied with the round discharge tip of the discharge electrode approaching the EUT as fast as possible (without causing mechanical damage) to touch the EUT. After each discharge, the ESD generator was removed from the EUT and re-triggered for a new single discharge. The test was repeated until all discharges were complete.

Vertical Coupling Plane (VCP):

The coupling plane, of dimensions 0.5m x 0.5m, is placed parallel to, and positioned at a distance 0.1m from, the EUT, with the Discharge Electrode touching the coupling plane.

The four faces of the EUT will be performed with electrostatic discharge.

Horizontal Coupling Plane (HCP):

The coupling plane is placed under to the EUT. The generator shall be positioned vertically at a distance of 0.1m from the EUT, with the Discharge Electrode touching the coupling plane.

The four faces of the EUT will be performed with electrostatic discharge.

2. Air discharges at insulation surfaces of the EUT.

It was at least ten single discharges with positive and negative at the same selected point.

5.4 TEST RESULT

Temperature:	22°C	Relative Humidity:	48%
Test Voltage:	DC 3.0V	Pressure:	1010hPa
Test Mode:	Running		


Mode	Air Discharge								Contact Discharge								Perform Criteria	Result
	4		8		10		15		2		4		6		8			
Test Location	+	-	+	-	+	-	+	-	+	-	+	-	+	-	+	-		
HCP									A	A	A	A					B	PASS
VCP									A	A	A	A						PASS
Slots	A	A	A	A														PASS
Surface	A	A	A	A														PASS

## 6 RADIATED, RADIO-FREQUENCY, ELECTROMAGNETIC FIELD IMMUNITY TEST (RS)

### 6.1 TEST SPECIFICATION

Basic Standard:	IEC/EN 61000-4-3
Required Performance:	A
Frequency Range:	80 MHz to 1000 MHz, 1400 MHz to 2000 MHz, 2000 MHz to 2700 MHz
Field Strength:	80 MHz to 1000 MHz: 3 V/m 1400 MHz to 2000 MHz: 3 V/m 2000 MHz to 2700 MHz: 1 V/m
Modulation:	1kHz Sine Wave, 80%, AM Modulation
Frequency Step:	1 % of fundamental
Polarity of Antenna:	Horizontal and Vertical
Test Distance:	3 m
Antenna Height:	1.5 m
Dwell Time:	$1.5 \times 10^{-3}$ decade/s

### 6.2 TEST SETUP

Note:

#### TABLE-TOP EQUIPMENT

The EUT installed in a representative system as described in section 7 of IEC/EN 61000-4-3 was placed on a non-conductive table 0.8 meters in height. The system under test was connected to the power and signal wire according to relevant installation instructions.

#### FLOOR-STANDING EQUIPMENT

The EUT installed in a representative system as described in section 7 of IEC/EN 61000-4-3 was placed on a non-conductive wood support 0.1 meters in height. The system under test was connected to the power and signal wire according to relevant installation instructions.

### 6.3 TEST PROCEDURE

The EUT and support equipment, which are placed on a table that is 0.8 meter above ground and the testing was performed in a fully-anechoic chamber.

The testing distance from antenna to the EUT was 3 meters.

The other condition need as following manners:

1. The frequency range is swept from 80 MHz to 1000 MHz, 1800 MHz, 2600 MHz, 3500 MHz, 5000 MHz, with the signal 80% amplitude modulated with a 1kHz sine wave. The rate of sweep did not exceed  $1.5 \times 10^{-3}$  decade/s. Where the frequency range is swept incrementally, the step size was 1% of fundamental.
2. The dwell time at each frequency shall be not less than the time necessary for the EUT to be able to respond.
3. The test was performed with the EUT exposed to both vertically and horizontally polarized fields on each of the four sides.

6.4 TEST RESULT

Temperature:	22°C	Relative Humidity:	48%
Test Voltage:	DC 3.0V	Pressure:	1010hPa
Test Mode:	Running		






Frequency Range (MHz)	RF Field Position	R.F. Field Strength	Azimuth	Perform Criteria	Result
80~1000	H / V	3 V/m (rms) AM Modulated 1000Hz, 80%	Front	A	PASS
			Rear		
			Left		
			Right		
1400~2000	H / V	3 V/m (rms) AM Modulated 1000Hz, 80%	Front	A	PASS
			Rear		
			Left		
			Right		
2000~2700	H / V	1 V/m (rms) AM Modulated 1000Hz, 80%	Front	A	PASS
			Rear		
			Left		
			Right		

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